propex•pixi™

apex locator

+
**WE
KNOW
ENDO.**

User Manual
A1030 000 000 00

CONTENTS

GB	English	5
F	Français	25
E	Español	47
P	Português	69
RU	Русский	91
CN	中国	117

Table of contents

Introduction	6
1. Indications for use	6
2. Contraindications	7
3. Warnings	7
4. Precautions	8
5. Adverse reactions	9
6. Step-by-step instructions	10
6.1 Content	10
6.2 Connecting the AC plug adapter	11
6.3 Recharging the battery	11
6.4 Replacement of the rechargeable battery	12
6.5 Cable connection test	13
6.6 Apex localization	14
6.7 Sound adjustment	16
6.8 Demo mode	16
6.9 Automatic shutdown	17
7. Maintenance of your propex-pixi™	17
7.1 General recommendations	17
7.2 Disinfection and sterilization procedure for lip clip, connection hook and fork	18
8. Troubleshooting	21
9. Warranty	23
10. Disposal of the product	24
11. Certification	24
12. Technical specifications	24
13. Standard symbols	24
APPENDIX	
Electromagnetic compatibility (EMC)	137

DIRECTIONS FOR USE

INTRODUCTION

Congratulations on the purchase of **propex-pixi™** apex locator. **propex-pixi™** is a device aimed at detecting the minor apical foramen based on analysis of electrical properties of different tissues inside the root canal system. For optimal safety and performance, read this user manual carefully before use. Make sure you have understood and followed the clinical precautions - as well as the general warnings, precautions and contraindications - before proceeding to determining a working length. Keep this user manual for future reference.

1. INDICATIONS FOR USE

propex-pixi™ is an electronic device used for apex location during root canal treatment. **propex-pixi™** must only be used in hospital environments, clinics or dental offices, by qualified practitioners.

2. CONTRAINDICATIONS

propex.pixi™ is not recommended for use:

- in patients who have a pacemaker or other implanted electrical devices, or have been cautioned by their physicians against the use of small electric appliances such as shavers, hair dryers, etc,
- in patients allergic to metals,
- in children.

3. WARNINGS

- The scale indication on the **propex.pixi™** screen does not represent a distinct length or distance in mm or other linear units. It simply indicates the file progression towards the apex.
- The following patient related factors may prevent accurate readings:
 - Blocked root canals,
 - Teeth with large apices,
 - Root fracture or perforation,
 - Metal crowns or bridges, if they come into contact with the file or the lip clip.
- Inaccurate or incorrect readings due to the environment are likely to occur in the following cases:
 - Presence of portable or movable radio frequency transmitters in the surroundings.
 - Film viewers or other Illumination devices which use an inverter may cause abnormal operation of the apex locator. Such devices should be turned off during use of the **propex.pixi™**.
 - Electromagnetic interference could cause improper operation of the device. In such cases the device behaviour may become abnormal or random. Usage of any devices emitting electromagnetic radiation, such as cellular phones, remote controls, transceivers, etc., should be prohibited in the vicinity of **propex.pixi™**.
- General safety warnings:
 - In order to prevent infectious agent transfer it is highly recommended to use a rubber dam system during the endodontic procedure.
 - Make sure that the lip clip, hook or fork does not come into contact with an electric power source such as an electrical socket. This could result in a severe electrical shock.
 - Do not use **propex.pixi™** in the presence of flammable substances.
- Only use the original charger.
- Use AAA 1.2V 1000mAh NiMH rechargeable batteries only. Usage of non-rechargeable batteries may cause device damage.

Recommended battery type: Manufacturer – GP Batteries
Model No.: GP100AAAHC.

4. PRECAUTIONS

Important notice:

The use of apex locators alone without a preoperative and postoperative radiograph is not recommended practice, since apex locators may not be able to work properly in all conditions. It is recommended to take an X-ray prior to the use of the unit, and compare the information given by both means. In case of clinically relevant discrepancy, it is recommended to take a second radiograph with a file in place at the selected working length [European Society of Endodontology. Quality guidelines for endodontic treatment: consensus report of the European Society of Endodontology. International Endodontic Journal, 39, 921–930, 2006].

It is important to follow the precautions below and pay close attention to any condition or situation that may influence the electrical conductivity during the procedure.

- Inaccurate or incorrect readings are likely to occur in the following cases, all procedure related:
 - Partially blocked canal.
 - Size of the measuring file differing significantly from the canal diameter. Ideally, the selected file should be the thickest one capable of reaching the apex.
 - Presence of liquids and/or tissue debris in the access cavity. Prior to the use of the device, the access cavity must be dried with a cotton pellet in order to prevent leaking current.
 - Contact of the file or the lip clip with metallic dental structures. Be particularly careful with patients fitted with metal crowns or bridges.
 - Contact of the file with another instrument.
 - Very dry canal, for instance in the presence of restoration. In this case the canal must be moistened with an irrigation solution, or with Glyde™ file prep.
 - Contact between the file and the gums (this may cause a false reading indicating that the apex has been reached).
 - Use of an ultrasonic scaler with the counter electrode attached to the patient (electrical noise from the scaler could interfere with the apex localization).
 - Use of the apex locator in conjunction with an electric scalpel.
 - Use of a damaged lip clip, hook or fork.
- For apex localization, concentrations of NaOCl higher than 5 % may result in reduced accuracy.
- As a safety precaution in order to avoid over-instrumentation, it is recommended to proceed as follows: place the file onto an endodontic ruler, where the **propex.pixi™** indicates '0.0'. Subtract a minimum of 0.5 mm from the measured file length.
- Please also respect the following precautions:
 - For your own safety, be aware of wearing personal protective equipment (gloves, glasses, mask).
 - If the bar graph makes sudden large movements in the coronal part of the canal, slowly continue advancing the file toward the apex until the signal returns to normal.

- This **propex.pixi™** unit must not be connected to or used in combination with any other apparatus or system. It must not be used as an integral component of any other apparatus or system. Using replacement parts or accessories not supplied by the original manufacturer or vender could adversely affect the EMC performance of the **propex.pixi™**.
- The device must be used with the manufacturer's original accessories only.
- Unplug the device before replacing the battery.
- Never use batteries that are leaky, deformed, discolored or otherwise abnormal.
- In case of battery leakage, carefully dry the battery terminals and remove all of the leaked liquid. Then replace the battery with a new one.
- Dispose of old batteries according to local codes and regulations.
- Accessories including lip clips, hooks or forks should be clean and without residue of chemical disinfectants or other medicinal solutions such as sodium hypochlorite or formalin.
- Do not expose **propex.pixi™** to any liquid.
- **propex.pixi™** must be stored in normal temperature (< 60°C) and humidity conditions.

5. ADVERSE REACTIONS

If the apex locator provides an incorrect reading and there is no radiographic data (see "Important notice" in the "Precautions" section in chapter 4), the following adverse reactions may occur:

- Incomplete root canal treatment,
- Apex perforation.

6. STEP-BY-STEP INSTRUCTIONS

6.1 Content

Check the content of the equipment before use:

1. **propex.pixi™** apex locator
 2. switching charger
 3. measurement cable with clip
 4. lip clip (2x)
 5. connection hook
- user manual
 - product card

6.2 Connecting the AC plug adapter

Select the plug adapter that matches your electric power outlet.

Pic. 1 Plug adapters for power supply

Align and insert the plug adapter at the rounded edge and then snap it into place by inserting the opposite end (see Pic. 1). To remove, pull the locking button (A) and pull out the plug adapter.”

6.3 Recharging the battery

propex.pixi™ is equipped with a rechargeable battery.

When the battery is low, a battery indicator appears on the status bar of the device. When the battery icon is blinking, the battery requires recharging. However, it is still functional for several treatments before the device shuts down.

Procedure for recharging the battery:

1. Disconnect the measurement cable
2. Connect the charger cable to **propex.pixi™**
3. Connect the charger to the mains

During battery charging the battery symbol will first blink (see Pic.2), then will remain steady when charging is completed (see Pic.3.)

Pic. 2 Charging

Pic. 3 Full

Duration of charging: Approximately 12 hours (24 hours after long period of non-use).

Note: **propex·pixi™** cannot be used while charging.

6.4 Replacement of the rechargeable battery

propex·pixi™ is powered by one 1.2V AAA NiMH rechargeable battery.

If a fully charged battery is not sufficient for normal device operation during at least one working day, the battery should be replaced by a new one at your earliest convenience. Please note that a new battery should be charged for 24 hours prior to the first use of the device.

Warning: Use GP100AAAHC or compatible NiMH rechargeable battery only. Usage of non-rechargeable batteries may cause device damage.

The battery compartment is located at the rear side of **propex·pixi™**.

1. Carefully lift and tilt back the silicon cover to expose the screw. Release the screw.

Pic. 4

2. Remove the battery compartment cover and the old battery.

Pic. 5

3. Insert a new battery into the battery compartment according to polarity marking.
4. Close the battery compartment, secure the cover with the screw and adjust the silicon screw cover.

The new battery should be charged for 24 hours prior to the first use of the device.

6.5 Cable connection test

A connection test feature is included in **propex.pixi™** in order to check the cables:

1. Connect the measurement cable and turn on the device.
2. Connect the metal part of the connection hook to the lip clip. Make sure that the accessories are cleaned properly before the test.
3. "Connection test" icon ➡⚡⚡⚡⚡⚡➡ should appear on the status bar - see Pic. 6.
4. If no icon appears, the connection hook or the measurement cable should be replaced.

Pic. 6

6.6 Apex localization

6.6.1 Getting started

Disconnect the charger from the device if connected.

1. Connect the measurement cable and turn on the device by pressing the ① “ON / OFF” button on the top of the device. The first bar will start blinking.
2. Attach the lip clip to the patient.
3. Insert the file into the canal.

Note: to ensure optimal performance the file size should be adjusted to the canal diameter.

4. Connect the connection hook to the metal shaft of the file.

The first bar will stop blinking accompanied by a double beep signal – see Pic. 7.

Pic. 7

Note: A blinking bar indicates faulty connection. Should this occur, check cable connections, clean the connection hook and the lip clip, moisten the canal if necessary and start again.

No other adjustments are required before starting apex localization.

6.6.2 Apex localization

Advance the file with slow clockwise turns. In the pre-apical zone bar 2.0 turns on – see Pic. 8 and an audio signal sounds. As the file progresses in the canal, subsequent bars turn on gradually (Pic. 9) and the interval between the audio signals becomes shorter. If the bar graph suddenly makes a large movement in the upper part of the canal, continue slightly towards the apex so the signal returns to normal.

Warning: The scale indication on the propex-pixi™ screen does not represent a distinct length or distance in mm or other linear units. It simply indicates the file progression towards the apex.

Pic. 8

Pic. 9

6.6.3 Apical zone

The apical zone is divided into 3 bars graduated from 1.0 to 0.0 (apex) – see Pic. 10-12.

Beginning of apical zone

Pic. 10

Middle of apical zone

Pic. 11

Apical position

Pic. 12

When the apex is reached, a solid tone is emitted.

The indication 0.0 on the **propex-pixi™** screen relates to the minor apical foramen file position (the apical length).

Note: As a safety precaution in order to avoid over-instrumentation, it is recommended to proceed as follows: place the file onto an endodontic ruler where the **propex-pixi™** indicates '0.0'. Subtract a minimum of 0.5 mm from the measured file length.

6.6.4 Over-instrumentation


A red “OVER” segment and an audio warning signal (rapid intermittent signal) indicate that the file has passed the apex – pic. 13.


Pic. 13

6.7 Sound adjustment

propex.pixi™ is equipped with an audio indicator which enables monitoring of the progression of the file within the canal in addition to visual monitoring.

The volume can be adjusted to four different levels: mute, low, normal and high, by successively pressing the  “VOLUME” button.

When sound level is muted the  icon is turned off. At other sound levels the  icon remains turned on.

6.8 Demo mode

The built-in demo mode is available to become acquainted with the device and to demonstrate its operation.

1. Disconnect the measuring cable or the charger from the device, if connected, and turn the device off.
2. To start the demo mode, press and hold the  “ON / OFF” button for about 2 seconds. until the whole display turns on (Pic. 14) and immediately turns off again accompanied by two beeps.

Pic. 14

3. During the demo cycle the device operating sequence is shown on the screen.
4. Demo cycles are repeated automatically until interrupted by the operator.
5. To exit the demo mode press the ⓘ "ON / OFF" button and hold for about 1 second until a beep sounds. The device shuts down.

Note: If the measurement cable is connected to **propex•pixi™** during the demo cycle, the device switches automatically to normal operation mode.

6.9 Automatic shutdown

propex•pixi™ automatically shuts down after 3 minutes of non-use. In order to prolong battery life, it is recommended to switch the device off after usage by pressing the ⓘ "ON / OFF" button.

7. MAINTENANCE OF YOUR PROPEX•PIXI™

7.1 General recommendations

- The device does not contain user serviceable parts. The service and repair should be provided by factory trained service personnel only.
- After each use, all objects that were in contact with infectious agents should be cleaned using towels impregnated with a disinfecting and detergent solution (a bactericidal, fungicidal and aldehyde free solution). Use of chemical agents may cause damage to the equipment. We recommend to use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval).

- The lip clip and the hook must be sterilized between treatments. Please note that the measuring cable cannot be autoclaved.
- In addition, a fork is not included in the **propex-pixi™** but it can be used and should follow the same procedure as the lip clip and hook.
- Follow the “Disinfection and sterilization procedure” described in section 7.2.
- The user is responsible for the sterility of the lip clip, the connection hook and the fork for the first cycle and each further usage.
- All damaged accessories should be discarded and dirty accessories should be cleaned and sterilized per the procedure described in section 7.2

7.2 Disinfection and sterilization procedure for lip clip, connection hook and fork

Foreword

For hygiene and sanitary safety purposes, the lip clip, the connection hook and the fork must be cleaned, disinfected and sterilized before each usage to prevent any cross-contamination between patients. This concerns the first use as well as the subsequent uses.

General recommendations

1 - Use only a disinfecting solution which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the disinfecting solution manufacturer. For all metal instruments, it is recommended to use anticorrosion disinfecting and cleaning agents

2 - For your own safety, please wear personal protective equipment (gloves, glasses, mask).

3 - The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments where applicable after sterility.

4 - Limitations and restrictions on reprocessing:

the appearance of defects such as cracks, deformations (bent, twisted), corrosion, loss of colour coding or marking, are indications that the devices are not able to fulfil the intended use with the required safety level.

5- Use only clean water in all cleaning and rinsing steps.

Step-by-step procedure

	Operation	Operating mode	Warning
1.	Disassembling	- Disassemble the device	
2.	Pre-Disinfection	- Soak all instruments immediately after use in a detergent and disinfecting solution combined with proteolytic enzyme if possible.	<ul style="list-style-type: none"> - Follow instructions and observe concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments). - The disinfecting solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as a corrosion inhibitor. - Do not use disinfecting solutions containing phenol or any products which are not compatible with the instruments (See general recommendations). - For visible impurities observed on instruments a pre-cleaning is recommended, brush manually with soft material.
3.	Rinsing	- Abundant rinsing (at least 1 min)	<ul style="list-style-type: none"> - Use clean water. - If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the cleaning.
4a.	Automated Cleaning with washer-disinfector	<ul style="list-style-type: none"> - Place the devices in a support or container to avoid any contact between instruments or posts. - Put them in the washer-disinfector (Ao value > 3000 or, at least 5 min at 90 °C). 	<ul style="list-style-type: none"> - Discard any instruments with large obvious defects (broken, bent). - Avoid any contact between instruments or posts when placing in the washer-disinfector. Use supports or container. - Follow instructions and observe concentrations given by the manufacturer (see also general recommendations). - Use only approved washer-disinfector according to EN ISO 15883, maintain and calibrate it regularly.
OR			
4b	Manual Cleaning or assisted by an ultrasonic device	<ul style="list-style-type: none"> - Place the devices in a support or container to avoid any contact between instruments. - Immerse in the disinfecting solution with cleaning properties, assisted by an ultrasonic device if suitable. 	<ul style="list-style-type: none"> - No visible impurities should be observed on the instruments. - Discard any instruments with large obvious defects (broken, bent, and twisted). - Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendations). - The disinfecting solution should be aldehyde free and without di- or triethanolamines as a corrosion inhibitor.

	Operation	Operating mode	Warning
5.	Rinsing	- Abundant rinsing (at least 1 min)	<ul style="list-style-type: none"> - Use quality water in accordance with local regulations. - If a disinfecting solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before the autoclaving. - Dry on a single use non-woven cloth, or with a drying machine or filtered compressed air.
6.	Inspection	<ul style="list-style-type: none"> - Inspect devices and sort out those with defects. - Assemble the devices (stops) 	<ul style="list-style-type: none"> - Dirty instruments must be cleaned and disinfected again. - Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument
7.	Packaging	- Place the devices in a support or container to avoid any contact between instruments or posts and pack the devices in "Sterilization pouches".	<ul style="list-style-type: none"> - Avoid any contact between instruments during sterilization. Use supports or containers. - Check the validity period of the pouch given by the manufacturer to determine the shelf life. - Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.
8.	Sterilization	- Steam sterilization at: 134 °C / 273°F, 18 min.	<ul style="list-style-type: none"> - The instruments, posts and the plastic supports must be sterilized according to the packaging labelling. - Use only autoclaves that match the requirements of EN 13060, EN 285. - Respect the maintenance procedure of the autoclave device given by the manufacturer. - Use only this recommended sterilization procedure. - Control the efficiency (packaging integrity, no humidity, colour change of sterilization indicators, physico-chemical integrators, digital records of cycles parameters).
9.	Storage	- Keep devices in sterilization packaging in a dry and clean environment	<ul style="list-style-type: none"> - Sterility cannot be guaranteed if packaging is open, damaged or wet. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and validity period).

8. TROUBLESHOOTING

Please review the checklist below should you experience a problem with your **propex-pixi™**. If the problem persists after following the proposed solutions, please contact your distributor.

Warning: The following patient related factors may prevent accurate readings:

- Blocked root canals,
- Teeth with large apices,
- Root fracture or perforation,
- Metal crowns or bridges, if they come into contact with the file or the lip clip.

#	Problem	Possible cause	Solution
1	<p>During battery charging the battery symbol blinks at a fast rate.</p> 	<p>1. The battery is not connected.</p> <p>2. The battery is not rechargeable.</p>	<p>1. Open the battery compartment, and connect the battery as described in the User Manual, section 6.4.</p> <p>2. Replace the battery with a rechargeable type as described in the user manual, section 6.4.</p>
2	<p>The device does not turn on by pressing the "ON / OFF" button.</p>	<p>1. The button is functioning incorrectly.</p> <p>2. The battery is discharged.</p> <p>3. Electronic malfunction.</p>	<p>1. Try pressing the "ON / OFF" button several times. If the device still does not turn on, please contact your distributor.</p> <p>2. Charge the battery.</p> <p>3. Contact your distributor.</p>
3	<p>The device shuts off during the procedure.</p>	<p>1. The battery is low.</p>	<p>1. Charge the battery.</p>
4	<p>No sound during the procedure.</p>	<p>1. The sound control is set at "Mute" level.</p>	<p>1. Adjust the sound level by pressing the "VOLUME" button.</p>

#	Problem	Possible cause	Solution
5	Display is not steady during the procedure.	<ol style="list-style-type: none"> 1. There is not a good contact between the lip clip and the oral mucosa. 2. The connection hook is soiled. 3. Deep caries provides a conductive path outside the canal. 4. Perforation. 5. Large lateral canal. 	<ol style="list-style-type: none"> 1. Ensure a good contact between the mucosa and lip clip (Place the lip clip in the labial angle opposite the tooth to be treated). 2. Clean the connection hook (with Ethanol). 3. Block the external conductive path. 4. Remove the file, close the perforation and repeat the apex detection procedure, carefully inserting the file into canal. 5. Try continuing the procedure by gently advancing the file.
6	The transmission of the electric signal is interrupted. The device does not show file progression inside the canal.	<ol style="list-style-type: none"> 1. Bad electrical contact. 2. The connection hook is not properly connected to the file. 3. The root canal is obliterated. 4. In the case of re-treatment: old filling material residues may block the root canal. 5. The root canal may be blocked by the remnants of a medication (e.g. calcium hydroxide). 	<ol style="list-style-type: none"> 1. Perform the cable connection test as described in the user manual, section 6.5. 2. Place the connection hook on the metal part of the file below the plastic handle. 3. Check the comparative X-ray image for hints. 4. Remove old root filling material residues prior to use. 5. Completely remove the remnants prior to use.

#	Problem	Possible cause	Solution
		<p>6. Root canal is extremely dry.</p> <p>7. The selected file is too small for a large root canal.</p> <p>8. Electronic malfunction.</p>	<p>6. Rinse the root canal with NaCl solution. Dry the access cavity with a cotton pellet/ air-blower.</p> <p>7. If there is no parietal contact use larger ISO size file. <u>Important:</u> exactly fitting files lead to precise results.</p> <p>8. Contact your distributor.</p>
7	Display reaction is erratic : “0.0” or “OVER” appear on the screen before the apical constriction is reached.	<p>1. Short circuit due to excess liquid (irrigation solution, saliva, blood) in the pulp chamber.</p> <p>2. A direct contact of the file with the gingiva or gingival proliferations, e.g. a fractured metal crown.</p> <p>3. A direct contact of the file with metal restorations (crown, parapulpal post, amalgam filling).</p>	<p>1. Dry the access cavity with a cotton pellet / air-blower. In case of excess bleeding wait until it has stops.</p> <p>2. For isolation: - adequate preparation filling of access cavity. - use a rubber dam.</p> <p>3. Isolate the file by inserting it into a small polyvinyl tube before use.</p>

9. WARRANTY

propex-pixi™ is warranted for 12 months from the date of purchase. The accessories (cables, battery etc.) are warranted for 6 months from the date of purchase.

The warranty is valid for normal usage conditions. Any modification or accidental damage will void the warranty.

10. DISPOSAL OF THE PRODUCT

Recycling : PLEASE DO NOT THROW AWAY! This product and all its components must be recycled through your supplier.

11. CERTIFICATION

propex-pixi™ complies with IEC 60601-1 Safety and IEC 60601-1-2 EMC (Electromagnetic compatibility) standards.

12. TECHNICAL SPECIFICATIONS

propex-pixi™ electronic apex locator belongs to the following category of medical devices:

- Internally powered equipment (AAA 1.2V 1000mAh NiMH rechargeable battery)
- Type BF applied parts
- Not suitable for use in the presence of flammable anesthetic mixtures with air, oxygen or nitrous oxide
- Continuous operation
- Ingress of liquids – not protected
- Environmental conditions during transportation: temperature: -20°C to +60°C (0 to 140°F); relative humidity: 10% to 90%, non-condensing

Specifications:

Dimensions: 66 x 55 x 18 mm
Weight: 55 gr.
Type of screen: Color LEDs
Supply: AAA 1.2V 1000mAh NiMH rechargeable battery
Switching charger: Input: 100-240 V AC ~ 50-60 Hz
Output: 6V DC \pm 5%, 1000 mA

13. STANDARD SYMBOLS

On the device label appear standard symbols as follows:

Class II equipment

Type BF applied part

Follow instructions for use

Manufacturer

Direct current

Serial no.

Recycling : PLEASE DO NOT THROW AWAY! This product and all its components must be recycled through your supplier.

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ELECTROMAGNETIC COMPATIBILITY (EMC)

Changes or modifications to this product not expressly approved by the manufacturer may result in increased emissions or decreased immunity performance of the product and could cause EMC issues with this or other equipment. This product is designed and tested to comply with applicable regulations regarding EMC and shall be installed and put into service according to the EMC information stated as follows.

WARNING

Use of portable phones or other radio frequency (RF) emitting equipment near the product may cause unexpected or adverse operation.

WARNING

The product must not be used adjacent to, or stacked with, other equipment. If adjacent or stacked use is necessary, the product must be tested to verify normal operation in the configuration in which it is being used.

Compliant Cables and Accessories

WARNING

The use of accessories, transducers and cables other than those specified may result in increased emissions or decreased immunity performance of the product. The table below lists cables, transducers, and other applicable accessories for which the manufacturer claims EMC compliance.

NOTE: Any supplied accessories that do not affect EMC compliance are not listed.

Description	Length. max
Measuring Cable	1.3 m
Charger Cable	1.5 m

Guidance and manufacturer's declaration - electromagnetic emissions

The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The Product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration - electromagnetic immunity

The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.

Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity for not life support equipment

The Product is intended for use in the electromagnetic environment specified below. The customer or the user of the Product should assure that it is used in such an environment.

Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the Product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m 80 MHz to 2.5 GHz	$d = 0,35 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 0,7 \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$

Immunity tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
			<p>Where P is the maximum output power rating in the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency range b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<p>a Fixed strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To access the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Product is used exceeds the applicable RF compliance level above, the Product should b observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Product.</p>			
b over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the not life support equipment

The Product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Product as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 0.35 \sqrt{P}$	80 MHz to 800 MHz $d = 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d = 0.7 \sqrt{P}$
0,01	0,04	0,04	0,07
0,1	0,11	0,11	0,22
1	0,35	0,35	0,7
10	1,1	1,1	2,2
100	3,5	3,5	7

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the higher frequency range applies.

Note 1: At 80 MHz and 800 MHz, the separation distance fort the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

www.dentsplymaillefer.com

Dentsply Maillefer
Chemin du Verger 3
CH – 1338 Ballaigues
Switzerland